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- Hall CB, Weinberg GA, Iwane MK et al. The burden of respiratory syncytial virus infection in young children. N. Engl. J. Med.360(6),588–598 (2009).
- Nair H, Nokes DJ, Gessner BD et al. Global burden of acute lower respiratory infections due to respiratory syncytial virus in young children: a systematic review and meta-analysis. Lancet. 375(9725),1545–1555 (2010).
- Englund. J.A., P. A. Peidra, A. Jewell, K. Patel, B. B. Baxter and E. Whimbey. 1996. Rapid diagnosis of respiratory syncytial virus infections in immunocompromised adults. J. Clin. Microbiol. 34:1649-1653.
- Crowcroft, N.S., F. Cutts and M.C. Zambon. 1999. Respiratory syncytial virus: an underestimated cause of respiratory infection, with prospects for a vaccine. Commun Dis Public Health. 2: 234-241.
- Barenfanger, J., C. Drake, N. Leon, T. Mueller and T. Troutt. 2000. Clinical and financial benefits of rapid detection of respiratory viruses: an outcomes study. J. Clin. Microbiol. 38:2824-2828.
- Biosafety in Microbiological and Biomedical Laboratories, 5th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (2007).
- Christopher King and Fred M. Henretig, Textbook of Pediatric Emergency Procedures, Chapter 122: 1125-1146, Obtaining Biologic Specimens. (1997).

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Vstrip RSV Rapid Test is for the detection of respiratory syncytial virus (RSV) antigen in patient's specimen by using nasopharyngeal (NP) swab and nasopharyngeal (NP) wash/ aspirate.

INTENDED USE

The Vstrip RSV Rapid Test is an immunochromatographic test for rapid, qualitative detection of RSV from NP swab and NP wash / aspirate specimen. This test is intended for in vitro diagnostic use to aid in the diagnosis of RSV infections in patients. It is recommended that negative test results be confirmed by cell culture or DFA.

SUMMARY AND EXPLANATION

Respiratory syncytial virus (RSV) is the most frequent cause of infantile bronchiolitis and pneumonia.^{1,2} Initial RSV infection often occurs in 6 weeks to 2 years old child, and the age of bronchiolitis and pneumonia caused by RSV infection is about 2 months to 6 months. Babies who are premature, congenital heart disease, bronchial lung hypoplasia, chronic lung disease, neuromuscular diseases, are considered as the high risk groups and requires more attention. RSV has also been implicated in severe respiratory infections in the elderly and immunocompromised.^{3,4} Since RSV infection may lead to serious or even fatal lung infection, therefore hospitalization is required for high-risk patients. Early symptoms of RSV infection may show symptoms such as fever and runny nose. When the virus spreads to the lungs, the symptoms will be aggravated which can be coughing, difficulty breathing, rapid breathing, wheezing. Since the early symptoms are similar to a common cold, it is hard to confirm the patient is RSV infected by these symptoms, therefore the use of Rapid Test can help physicians have early identification of the virus, give the correct dosage to patients, and prevent unnecessary antibiotics or other treatments being used.⁵ The Vstrip RSV Rapid Test is an immunochromatographic test that uses antibody-coated colloidal gold to detect the presence of RSV antigen in NP swab and NP wash / aspirate from symptomatic patients. The test is easy to perform and test results can be visually interpreted in 10 minutes.

PRINCIPLE OF THE TEST

The Vstrip RSV Rapid Test is a lateral flow immunochromatographic assay for the direct visual detection of RSV antigen in extracted specimen. When extracted specimen is added to the sample well, RSV antigen binds to the antibody coupled immunogold in the strip to form an antigen-antibody complex. This complex migrates across the test strip membrane and is captured by the line of RSV monoclonal antibody to form an antibody-antigen-antibody sandwich and a reddish-purple line is observed in the test region. An independent 2nd functional control line in the control region will always appear which indicates that the test has been performed correctly. If RSV antigen is not present, or is present at very low level in the specimen, only the control line will be visible. If the control line does not develop, the test is invalid.

MATERIALS AND REAGENTS PROVIDED

1. Test cassettes (20) – Mouse monoclonal anti-RSV antibodies.
2. Extraction buffer vial (20) – with detergent, protein and salt.
3. Disposable droppers (20)
4. Sterile nasopharyngeal swabs (20)
5. Positive control reagent (1): With non-infectious RSV antigen, detergents, protein and salt.
6. Negative control reagent (1): With detergents, protein and salt.
7. Package Insert (1)

MATERIALS NOT PROVIDED

1. Specimen collection container
2. Timer
3. Disposable gloves

WARNINGS AND PRECAUTIONS

1. For *in vitro* diagnostic use.
2. Directions should be read and followed carefully.
3. Do not use the kit contents beyond the expiration date printed on the outside of the box.
4. Do not interchange or mix different lots of reagents.
5. Do not add more than specified volume (3 drops) of sample extract buffer or positive/negative control to cassette.
6. Disregard test results beyond specified time (10 min).
7. Use appropriate precautions in the collection, handling, storage, and disposal of patient samples and used kit contents.⁶
8. Use of Nitrile or Latex gloves is recommended when handling patient samples.⁶
9. Dispose of containers and used contents in accordance with Federal, State and Local requirements.
10. Do not reuse kit components or test devices.
11. The Test Cassettes must remain sealed in the protective foil pouch until use.

12. To obtain accurate results, you must follow the Package Insert.
13. Inadequate or inappropriate specimen collection, storage, and transport may yield false test results.
14. Seek specific training or guidance if you are not experienced with specimen collection and handling procedures.

STORAGE INSTRUCTION

1. The product should be stored at 15-30°C, away from direct sunlight.
2. Do not freeze or overheat the test kit or kit reagents.
3. Kit contents are stable until the expiration date printed on the outer box.
4. The test cassettes must be kept in the foil pouch until use.

SPECIMEN COLLECTION AND HANDLING

Collect patient specimen by following the recommended guideline⁷ or your institution's protocol for obtaining NP swab and NP wash / aspirate specimen. For NP wash / aspirate specimen, use the minimal amount of saline that your procedure allows, as excess volume will decrease the concentration of antigen in the specimen.

ASSAY PROCEDURE

Performing the assay outside the time and temperature ranges provided may produce invalid results. Assays not performed within the established time and temperature ranges must be repeated.

Expiration date:

Check expiration on each individual test package or outer box before using. Do not use any test past the expiration date on the label.

NP Swab Test Procedure

1. Place the swab with sample into the reagent tube. Roll the swab at least five (5) times while pressing the head of swab against the tube bottom and the side of the reagent tube. Leave the swab in the reagent tube for one (1) minute.
2. Roll the swab head against the inside wall of reagent tube as you remove it. Dispose of the used swab according to the biohazard waste disposal protocol.
3. Use the disposable dropper to transfer 3 drops of sample extract buffer from the reagent tube to the sample well of the cassette.
4. Read result at 10 minutes. Some positive results may appear sooner. Do not read the result after 10 minutes.

NP Wash/Aspirate Test Procedure

1. Use the disposable dropper to draw 0.5 ml of the NP wash/aspirate sample and add to the extraction tube.
2. Screw the cap of the extraction tube and swirl it to mix. Wait 1 minute to allow the mixture to react.
3. Use the disposable dropper (the same dropper of item 1) to transfer 3 drops of sample extract buffer from the reagent tube to the sample well of the cassette.
4. Read result at 10 minutes. Some positive results may appear sooner. Do not read the result after 10 minutes.

INTERPRETATION OF RESULTS

Positive Result


At 10 minutes, a visible reddish-purple line appears in the window next to the Test Line (T), and a reddish-purple line next to the Control Line (C). This indicates RSV antigen was detectable in the specimen.

Negative Result

At 10 minutes, the appearance of only one reddish-purple line in the Control Line (C) and no reddish-purple line in the Test Line (T) indicates that the sample is negative for RSV antigen or the antigen level is below the detection limit.

Invalid

If at 10 minutes, a reddish-purple line next to the Control (C) is not visible or the background color interferes with interpretation of the test or control lines, the test is uninterpretable. An uninterpretable test must be repeated, a new specimen obtained and retested, or the specimen sent to the clinical laboratory for culture isolation.

QUALITY CONTROL

The Positive Controls should be assayed once upon receipt of the kit. Internal procedural controls are included in the test. A reddish-purple line appearing in the control region is an internal control. It confirms sufficient specimen volume and correct procedural technique. The Positive Controls, a visible reddish-purple line appears in the window next to the Test Line (T). The positive controls are intended to monitor for reagent failure, but will not ensure precision at the analytical assay cutoff.

QUALITY CONTROL TESTING PROCEDURE

Positive/ Negative Control Test Procedure

1. Add 3 drops of control solution to the sample well on the cassette directly.
2. Read the result at 10 minutes. Do not read result after 10 minutes.

LIMITATIONS OF THE PROCEDURE

1. The contents of this kit are to be used for the qualitative detection of RSV antigen from the NP swab and NP wash / aspirate specimen.
2. A negative test result may occur if the level of antigen in a sample is below the detection limit of the test.
3. Failure to follow the test procedure and interpretations of test results may adversely affect test performance and/or invalidate the test results.
4. Test results must be evaluated in conjunction with other clinical data available to the physician.
5. Negative test results do not rule out other potential non-RSV viral infections. Negative results should be confirmed by cell culture.
6. Positive test results do not rule out co-infections with other pathogens.
7. Monoclonal antibodies may fail to detect, or detect with less sensitivity, RSV viruses that have undergone minor amino acid changes in the target epitope region.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity (Limit of Detection)

The analytical sensitivity was evaluated using 2 RSV strains. Each strain was serially diluted in Vstrip extraction buffer. Strains were assayed in triplicate using the Vstrip RSV Rapid Test until no positive signal could be observed. Results are summarized below:

Subgroup	RSV strain	TCID ₅₀ / mL
A	A2 (ATCC VR-1540)	1.12 x 10 ³
B	18537 (ATCC VR-1580)	8.89 x 10 ³

TCID₅₀/mL = 50% tissue culture infectious dose

Cross reactivity study

The Vstrip RSV Rapid Test was evaluated with a total of 15 bacteria, 7 viruses. Bacteria were evaluated at 10⁷ CFU/mL. Viruses were evaluated at a concentration of at least 10⁵–10⁷ TCID₅₀/mL. None of the microorganisms tested in the following table gave a positive result in the Vstrip RSV Rapid Test.

Bacterial Panel

<i>Pseudomonas aeruginosa</i>	<i>Serratia marcescens</i>
<i>Staphylococcus aureus</i>	<i>Staphylococcus epidermidis</i>
<i>Streptococcus Group A</i>	<i>Streptococcus Group B</i>
<i>Streptococcus Group C</i>	<i>Streptococcus Group F</i>
<i>Streptococcus Group G</i>	<i>Streptococcus mutans</i>
<i>Streptococcus pneumoniae</i>	<i>Streptococcus sanguis</i>
<i>Salmonella typhi</i>	<i>Salmonella enteritidis</i>
<i>Vibrio cholerae</i>	

Viral Panel

Adenovirus type 7	Human parainfluenza virus type 2
Echovirus type 11	Herpes simplex virus type2
Influenza A virus H1N1	Influenza A virus H3N2
Influenza B virus	

Interfering Substances

Hemoglobin, common chemicals and a commercial nasal spray product were evaluated and did not interfere with the Vstrip RSV Rapid Test at the levels tested below:

Interference substances/ Concentration

Hemoglobin (2%)	Oxymetazoline HCl (10 mg/mL)
Aspirin/ Acetylsalicylic Acid (20 mg/mL)	Phenylephrine HCl (100 mg/mL)
Dextromethorphan (10 mg/mL)	Swinin nasal sprays (10%)
Diphenhydramine HCl (5 mg/mL)	

PACKAGING

Vstrip RSV Rapid Test.....20 Tests/Kit

ORDERING INFORMATION

Product No. : IG02020C

REFERENCES

1. Hall CB, Weinberg GA, Iwane MK et al. The burden of respiratory syncytial virus infection in young children. N. Engl. J. Med.360(6),588–598 (2009).
2. Nair H, Nokes DJ, Gessner BD et al. Global burden of acute lower respiratory infections due to respiratory syncytial virus in young children: a systematic review and meta-analysis. Lancet375(9725),1545–1555 (2010).
3. Englund. J.A., P. A. Peidra, A. Jewell, K. Patel, B. B. Baxter and E. Whimbey. 1996. Rapid diagnosis of respiratory syncytial virus infections in immunocompromised adults. J. Clin. Microbiol. 34:1649-1653.
4. Crowcroft, N.S., F. Cutts and M.C. Zambon. 1999. Respiratory syncytial virus: an underestimated cause of respiratory infection, with prospects for a vaccine. Commun Dis Public Health. 2: 234-241.
5. Barenfanger, J., C. Drake, N. Leon, T. Mueller and T. Troutt. 2000. Clinical and financial benefits of rapid detection of respiratory viruses: an outcomes study. J. Clin. Microbiol. 38:2824-2828.
6. Biosafety in Microbiological and Biomedical Laboratories, 5th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (2007).
7. Christopher King and Fred M. Henretig, Textbook of Pediatric Emergency Procedures, Chapter 122: 1125-1146, Obtaining Biologic Specimens. (1997).

SYMBOL LEGEND

Catalog Number	Batch Code	Contains sufficient for <n> tests
In vitro diagnostic medical device	Positive Control Negative Control	Do not reuse
Consult instructions for use	Use-by date	CE Marking
Temperature Limit	Manufacturer	Authorized Representative in the European Community

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Panion & BF Biotech Inc.

16F., No.3, Park St., Nangang Dist., Taipei City 115, Taiwan

Manufactured by: Panion & BF Biotech Inc. at Xizhi

Manufacturer Address:

6F.-3, No.306, Sec. 1, Datong Rd., Xizhi Dist., New Taipei City

22146 Taiwan

Technical Support: +886-2-2691-9895

Email: xizhi@pbf.com.tw

www.pbf.com.tw

www.vstriptechnology.com

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HiSS Diagnostics GmbH, Tullastraße 70