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- Wannamaker, L.E. (1979) Reviews of Inf. Dis. 1:967-973 .
- Manual of Clinical Microbiology, 6th Edition, ASM. (1995) pg. 301.
- Maxted, W.R.J. (1953) Clinical Pathology. 6:224-226.
- Levinson, M.L., and P.F. Frank. (1955) J. Bacteriology. 69:284-287.
- Youmans, G.P., Paterson, P.Y., and Sommers H. M. (1975) The Biological and Clinical Basis of Infectious Disease. pp. 172-185.

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Immunochromatographic Rapid Test for the detection of Strep A carbohydrate antigens in throat swab specimens. For professional in vitro diagnostic use only.

INTENDED USE

The PBF Vstrip Strep A Rapid Test is a rapid chromatographic immunoassay for the qualitative detection for Group A Streptococcal antigen from throat swabs. The test is used to aid in the diagnosis of Group A Streptococcal infection.

SUMMARY AND EXPLANATION

Group A Streptococcal is the most significant pathogen causing pharyngitis. In the case of Group A Streptococcus, antibiotic therapy is the treatment of choice. If left untreated, serious such as rheumatic fever may occur. Conventional methods for detecting Group A Streptococcus infection involve 24-48 hour culture of throat swab specimens followed by confirmation of beta-hemolytic colonies as Group A Streptococcus^{2,3,4}. Immunological confirmation of Group A Streptococcus is based on the antigenic characteristics of the group specific carbohydrate⁵. The PBF Vstrip Strep A Rapid Test is a rapid test to qualitatively detect the presence of Group A Streptococcus antigen in throat swab specimens, providing results within 5 minutes.

PRINCIPLE OF THE TEST

The PBF Vstrip Strep A Rapid Test is an immunochromatographic membrane assay to detect Group A Streptococcus antigen from throat swab. To perform the test, a throat swab specimen is collected. Antigen is extracted from the swab specimen with reagents A and B. The dipstick is added to the extracted sample. If Group A Streptococcus is present in the sample, it will form complexes with the anti-Group A streptococcus antibody conjugated color particles located at two separate locations on the test stick. The complex will then be bound by the anti-Group A Streptococcus capture antibody and visible red test line will appear to indicate a positive result. A blue control line will also appear to indicate the test is valid.

REAGENTS AND MATERIALS

- 20 Test Dipsticks
- 20 Test Tubes
- 20 Throat Swabs
- 10 ml Extraction Reagent A (2M Sodium Nitrite)
- 10 ml Extraction Reagent B (0.4M Acetic Acid)
- 1 mL Positive Control (Non-viable Strep A containing preservative)
- 1 mL Negative Control (Non-viable Strep C containing preservative)
- Package Insert

MATERIALS NOT SUPPLIED

- Timer or watch

WARNINGS AND PRECAUTIONS

1. For In Vitro Diagnostic Use.
2. Do not use beyond the expiration date
3. Dispose of containers with used contents in accordance to Federal, State and Local requirements.
4. The dipstick must remain sealed in the protective foil pouch until just prior to use.
5. Do not eat, drink or smoke in the area where the specimens and kits are handled.
6. Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout the procedure and follow the standard procedures for proper disposal of specimens.
7. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are assayed.
8. WARNING: Extraction Reagent A is harmful if swallowed or absorbed through skin. May cause eye irritation.
9. CAUTION: Extraction Reagent B may cause skin, eye and respiratory tract irritation.
10. The positive and negative controls contain sodium azide (NaN₃) as a preservative.
11. Do not interchange reagent bottle caps.
12. Do not interchange external control solution bottle caps.
13. Humidity and temperature can adversely affect results.

STORAGE

- Store Test Sticks and reagents tightly capped at 15°C- 30°C.
- Do not use Test Sticks or reagents after expiration data.

SPECIMEN COLLECTION AND HANDLING

Using the swab provided in the kit (or use rayon tipped swabs to collect throat specimens. Do not use calcium alginate, cotton tipped or wooden shaft swabs.), collect specimen by sampling material from the back of the throat and any white patches in the tonsillar area. Avoid contact with teeth, gums, and tongue or cheek surfaces. Process patient sample immediately following specimen collection. If a culture result is also desired, inoculate the culture plate with the swab prior to performing the PBF Vstrip Strep A Rapid Test by rolling it firmly over a small portion of one quadrant of the agar plate. Do not perform the PBF Vstrip

Strep A Rapid Test before inoculating the culture plate, as the extraction reagents render the bacteria on swab nonviable. When immediate testing is not possible, place the patient swab into a dry test tube for transport or storage. Throat swab specimens may be stored refrigerated at 2-8°C for up to 72 hours. The use of transport media is not recommended with the product.

QUALITY CONTROL

Built-in Control Features

The PBF Vstrip Strep A Rapid Test provides three levels of internal procedural controls with each test run. For daily quality control, PBF recommends documenting that these internal controls were checked for the first sample tested each day. The color of the Extraction Reagent changes from pink to clear as the reagents are mixed together. The color change is an internal extraction reagent control and is an indication that the reagents were mixed and functioning properly. The appearance of a blue Control Line is an internal control. The Dipstick must absorb the proper amount of sample and the Dipstick must be working properly for the blue Control Line to appear. Additionally, the appearance of the Control Line indicates that capillary flow occurred. A clear background is an internal background negative control.

External Quality Control Testing

External controls may also be used to demonstrate that the reagents and assay procedure perform properly. PBF recommends that positive and negative controls be run once for each untrained operator, once for each new shipment of kits-provided that each different lot received in the shipment is tested- and as deemed additionally necessary by your internal quality control procedures, and in accordance with local, state, and federal regulations or accreditation requirements. If the controls do not perform as expected, repeat the test or contact PBF Technical Service. Support before testing patient specimens.

ASSAY PROCEDURE

All clinical specimens must be at room temperature before beginning the assay.

1. Just before testing, add 5 drops Reagent A (pink) and 5 drops Reagents B to the Test tube (the solution should turn light yellow).
2. Immediately put the swab into the Tube. Vigorously mix the solution by rotating the swab forcefully against the side of the Tube at least ten (10) times. Best results are obtained when the specimen is vigorously extracted in the solution. Let stand for 1 minutes.
3. Express as much liquid as possible from the swab by squeezing the sides of the tube as the swab is withdrawn. Discard the swab.
4. Remove the Test Stick(s) from foil pouch. Place the Dipstick into the tube with the arrows of the Dipstick pointing down.
5. Read result at five (5) minutes. Some positive results may appear sooner.

Interpretation of Results

Positive result: Two distinct lines

One line should be in the control region (C) and another line should be in the test region (T). A positive result indicates that Strep A was detected in the sample. NOTE: The intensity of the red color in the test line region (T) will vary depending on the concentration of Strep A present in the sample. Therefore, any visible shade of red in test region (T) should be considered positive.

Negative result: One line in the control region (C)


NO apparent red appears in test region (T). A negative result indicates that Strep A is not present in the sample, or is present below the detectable level of test. The patient's sample should be cultured to confirm the absence of Strep A infection. If clinical symptoms are not consistent with results, obtain another sample for culture.

Invalid: Control line fails to appear.

Insufficient sample volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test dipstick. If the problem persists, discontinue using the kit immediately and contact PBF Technical Service.

LIMITATIONS

- The PBF Strap Strep A Rapid Test does not differentiate between viable and nonviable organisms. Patient that have recently been treated for Strep A or

similar infections may positive results for a period of time following effective treatment due to the presence of Group A strep antigen in nonviable organism.

- Pharyngitis can be caused by organisms other than Group A Streptococcus. Further diagnostic testing, including culture, should be performed if laboratory findings are inconsistent with clinical presentation.
- Excess blood or mucus on the swab specimen may interfere with test performance and may yield a false positive result. Avoid touching the tongue, cheeks, and teeth and any bleeding areas of the mouth with the swab when collecting specimens.
- The PBF Vstrip Strep A will not differentiate asymptomatic carrier of Group A Streptococcus from those exhibiting streptococcal infection.
- A negative result may be obtained if the amount of extracted antigen is below the sensitivity of the test. Culture confirmation is recommended for all PBF Vstrip Strep A Rapid Test negative test results.
- A single swab that is used both to inoculate a culture plate and to perform the rapid test may have reduced sensitivity in the PBF Vstrip Strep A Rapid Test.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity (Limit of Detection)

All of Group A Streptococcus strain from BCRC (Taiwan) were tested as positive results in the PBF Vstrip Strep A Rapid Test. Strep A organisms were quantitated using standard microbiological methods. Analytical sensitivity studies, performed by testing dilutions of culture of Group A Streptococcus, ATCC 19615, determined the PBF Vstrip Strep A Rapid Test limit of detection to be 4.0 x10⁴ CFU per test.

Analytical specificity (Cross-Reactivity):

To determine the analytical specificity of the PBF Vstrip Strep A Rapid Test, 34 commensal and pathogenic microorganisms (14 bacteria and 20 viruses) that may be present in the pharynx were tested. All of the following microorganisms were negative when tested at concentrations ranging from 10² to 10⁹ TCID₅₀/ml (viruses) and 10⁷ to 10⁸ organisms/ml (bacteria).

Virus	Bacteria
Influenza A virus H1N1 (swine) Taiwan 2009	Pseudomonas aeruginosa
Influenza A virus H1N1 (swine) Taiwan 2010	Serratia marcescens
Influenza A virus H1N1 (human) WSN 1933	Staphylococcus aureus
Influenza A virus H1N1 (human) PR 1934	Staphylococcus epidermidis
Influenza A virus H1N1 (swine) Taiwan 2014	Streptococcus Group B
Influenza A virus H1N1 (swine) Taiwan 2011	Streptococcus Group C
Influenza A virus H3N2 Taiwan 2011	Streptococcus Group F
Influenza A virus H3N2 Taiwan 2013	Streptococcus Group G
Influenza A virus H3N2 Taiwan 2014	Streptococcus mutans
Influenza A virus H3N2 Taiwan 2015	Streptococcus pneumoniae
Influenza B virus Taiwan 2014-1	Streptococcus sanguis
Influenza B virus Taiwan 2014-2	Salmonella typhi
Influenza B virus Taiwan 2014-3	Salmonella enteritidis
Adenovirus type 7	Vibrio cholerae
Adenovirus type 41	
Respiratory syncytial virus 2958-14	
Respiratory syncytial virus 2313-15	
Human parainfluenza virus Type 2	
Echovirus Type 11	
Herpes simplex virus Type 2	

Interfering substances:

The following substances, naturally present in throat specimens or that may be artificially introduced into the pharynx, were evaluated in the PBF Vstrip Strep A Rapid Test at the concentrations listed below and were found not to affect test performance.

Substance	Concentration
Ibuprofen	10 mg/ml
Aspirin/Acetylsalicylic Acid	20 mg/ml
Dextromethorphan	10 mg/ml
Diphenhydramine	5 mg/ml
Oxymetazoline HCl	10 mg/ml
Phenylephrine HCl	100mg/ml
Swinin nasal spray	10%

Hemoglobin	2%
Mucin	1%

Clinical correlation

TBD

REFERENCES

1. Wannamaker, L.E. (1979) Reviews of Inf. Dis. 1:967-973 .
2. Manual of Clinical Microbiology, 6th Edition, ASM. (1995) pg. 301.
3. Maxted, W.R.J. (1953) Clinical Pathology. 6:224-226.
4. Levinson, M.L., and P.F. Frank. (1955) J. Bacteriology. 69:284-287.
5. Youmans, G.P., Paterson, P.Y., and Sommers H. M. (1975) The Biological and Clinical Basis of Infectious Disease. pp. 172-185.

PACKAGING

PBF Vstrip Strep A.....20 Tests/Kit

ORDERING INFORMATION

Product No. : IG03020S

SYMBOL LEGEND

	Catalog Number
	For In Vitro diagnostic use
	Consult instructions for use
	Temperature Limitation
	Batch Code
	Positive Control
	Negative Control
	Use By (Expiration Date)
	Manufacturer
	Non-reusable
	CE Marking

PBF & BF Biotech Inc.

Manufactured by:
 Panion & BF Biotech Inc. Xizhi Factory
 Manufacturer Address:
 6F-3, No.306, Sec. 1, Datong Rd.,
 Xizhi Dist., New Taipei City 22146
 Taiwan (R.O.C.)
 Technical Support : +886-2-26919895
 Email : xizhi@pbf.com.tw